

# **Atomization of CMP Slurry and Its Characterization**

**Handol Lee, Seong Chan Kim, Shawn Chen and David Y. H. Pui**

**Particle Technology Laboratory  
Mechanical Engineering  
University of Minnesota**

# Introduction


- **Chemical mechanical planarization (CMP) is a material processing technique used to polish semiconductors.**
- **Abrasive particles in CMP slurry have a **size distribution** which directly affects critical metrics including **rate of removal and wafer defects**.**
- **Particle size analysis is therefore a key indicator of CMP slurry performance.**

## Objectives

- **Characterization of CMP slurry, e.g., slurry particle **size distribution****



# SEM Images of CMP



- **CMP slurry particles on a 400 nm rated Nuclepore filter**
- **Small particles ( $< 10$  nm) and large aggregates are shown in SEM images.**



# Size Measurement Methods



- **Dynamic light scattering (DLS)**
  - ✓ Measure the intensity of the scattered light fluctuates
  - ✓ Detection limit of particle size:  $> 1\text{nm}$
- **Nanoparticle tracking analysis (NTA)**
  - ✓ Consist of laser and CCD camera to visualize particle behaviors
  - ✓ Measure particle size from diffusive Brownian motion of liquid-borne particles
  - ✓ Detection limit of particle size:  $> 20 \sim 30 \text{ nm}$
- **Aerosolization (scanning mobility particle sizer, SMPS)**
  - ✓ Dispersion by electrospray or atomizer
  - ✓ Measure airborne particle size distribution by SMPS
  - ✓ Detection limit of particle size:  $> 1 \text{ nm}$



# DLS Measurement of CMP



**Brownian motion**






- **Suspended particles exhibit Brownian motion.**
- **Particle velocity and size distributions are measured by light scattered from moving particles.**
- **Drawback**
  - ✓ **Low resolution (high concentration required)**
  - ✓ **Poor detection capability for particles with different sizes (mixtures)**
  - ✓ **Small particles hidden by larger particles**






# DLS Measurement of CMP with Different Concentrations

1<sup>st</sup> run




2<sup>nd</sup> run




- **Unstable (not consistent) size distribution, e.g., CMP 100 ppm**
- **Size distribution varies a lot with CMP concentration.**



# Nanoparticle Tracking Analysis (NTA)



**Recording**



**Tracking**

- **Composed of ultramicroscope, laser illumination and CCD camera**
- **Capturing light scattered by NPs using CCD camera**
- **Tracking Brownian motion of NPs from frame to frame**
- **Calculating particle size through the Stokes-Einstein equation**
- **Different setting required depending on types of samples (camera level and detection threshold)**






# NTA Measurement of CMP

**Concentration: ppb range (0.005 ppm)**

**Much lower concentration required compared to DLS**



- **Small particles (< 10 nm) cannot be seen by NTA.**
- **Visualized particles – larger than 20~30 nm (detection limit)**
- **Small particles (shown in the SEM image) are hidden by larger particles or cannot be seen (< 10 nm).**



UNIVERSITY OF MINNESOTA



# Aerosolization Methods

## DISPERSION



**Atomizer**



**Electrospray**

## MEASUREMENT



**Scanning mobility particle sizer (SMPS)**


- **Atomizer and electrospray are used for generating airborne particles.**
- **SMPS measures the size distribution of airborne particles (single particle counting).**
- **However, when using aerosolization method, residues contributed from water impurities and surfactant can interfere with the main particles.**





UNIVERSITY OF MINNESOTA

# CMP Size Distribution Obtained by Atomizer - SMPS

Only ultrapure water



CMP slurry





- **Residue particle size distribution from impurities of ultrapure water cannot be negligible (large droplets from atomizer).**
- **Due to the large droplets containing multiple particles, agglomeration occurs when the liquid is evaporated.**



# Dispersion Methods

## Electrospray (ES)



- **ES generates very small monodisperse droplets ( $D_d \sim 200$  nm).**
- **Due to the generation of small droplets, it can generate singlet particles (one particle in one droplet).**









# Mixtures of Different Particle Sizes

**DLS**

**DLS vs ES-SMPS**

**ES-SMPS**



- **DLS cannot distinguish the two particle sizes.**
- **DLS output is dominated by the larger particle size, i.e., 40 nm.**
- **ES-SMPS is based on single particle counting, which enables to show completely distinguished 10 and 40 nm Au particles.**




UNIVERSITY OF MINNESOTA


# Mixtures of Different Particle Sizes

## NTA vs ES-SMPS

**NTA**




**ES-SMPS**




- **NTA cannot distinguish two particles.**
- **NTA is also based on intensity measurement.**
- **The NTA output is dominated by larger particles.**



# Zeta Potential of CMP Particles



## Capillary tube in electrospray




- **CMP slurry particles in ultrapure water have positive 55.7 mV zeta potential at pH 3.5 ~ 4.**
- **However, the silica surface inside capillary tube of ES has negative zeta potential, resulting in deposition of CMP particles during transport.**



# Conclusion and Future Work

- **Intensity based measurements (DLS, NTA) show biased signals to larger particles.**
- **Aerosolization method is a promising technique to characterize polydisperse (or mixtures of different sizes) particles.**



- **In high pH ( $> 9$ ) solution, CMP slurry particles tend to have negative zeta potentials.**
- **No deposition inside capillary tube will be expected when particles have negative zeta potentials.**



# THANK YOU

## Q/A