

Measurement of Particle Geometric Surface Area Using a Modified Weighted Sum Method with a Time Resolution of 1S

Lipeng Su, Qisheng Ou, Nanying (Leo) Cao, and David Y. H. Pui

Center for Filtration Research, University of Minnesota

Presented at 53rd CFR Review Meeting, Prior Lake, MN

April 25th, 2018




Outline


- ☐ Background and objective
- ☐ Methodology
- ☐ Experiments
- ☐ Summary




Significance of surface area




Catalytic activity



Health effect



Drug action



Gas absorption



Main methods for surface area: offline


Brunauer–Emmett–Teller (BET)



Character:

- Surface area including pores
- Gas adsorption, direct
- High detection limit
- Ex-situ
- Time consuming, bulky & costly

Electron microscopy (TEM and SEM)




Character:

- 2D projected area-morphological information
- Low accuracy
- Ex-situ
- Time consuming, bulky & costly



Main Methods for surface area: online

Scanning Mobility Particle Sizer



- In minutes
- Working fluid
- Bulky & costly

Nanoparticle Surface Area Monitor



- Only lung-deposited surface area
- Reasonable size and cost
- 1s resolution

Electrical Low Pressure Impactor (+)



- Inertia separation
- Wide size stage
- Low accuracy
- costly





Main Method for surface area: online 2

Based on [NSAM](#), *Cao and Pui et al., (2017a, 2017b) developed a Geometric Surface Area Monitor (GSAM) and an weighted-sum method.

- However, GSAM Needs additional impactor even to spherical particle;
- Weighted-sum method Needs ten or more seconds of response time

$$I_{c1} + aI_{c2} \propto d_{p1}^2 + d_{p2}^2 + \dots d_{pn}^2$$



So far, **NO instrument** measures the particle geometric surface area (GSA) with a time-resolution of 1s



Objective

A simple method measuring geometric surface area (GSA) concentration of aerosol nanoparticles in real time with **1s resolution**.

Features:

- GSA, 1s response time
- Cost effective and Transportation friendly
- no working fluid, no radiation source, high accuracy, and large range.



Outline

☐ Background and objective

☒ Methodology

☐ Experiments



☐ Summary



Methodology

Single Channel (Cao and Pui)

Two Parallel Channels
(This Work)



UNIVERSITY OF MINNESOTA

Outline

☐ Background and objective

☐ Methodology


☒ Experiments

☐ Summary



1. Calibration (electrometer) setup

setup





- **Monodisperse Singly charged particle**
- **Reference electrometer** : Using electrometer TSI 3068B



1. Calibration (electrometers)

results




- Testing particle size: 40nm, 100nm, 150nm, 200nm, 250nm, 300nm
- Good accuracy could be found in both electrometers
- Reference electrometer : Using electrometer TSI 3068B



2. Calibration (Sensitivity)

setup




- **Ion trap: low voltage**
- **High voltage: from 300V to 2400V**



2. Selection of appropriate voltage (spheres)



- **Low voltage: 20V (default voltage):** only removing excess ions; covering the whole size range
- **High voltage: Which voltage should we choose?**



$$S = I/C$$

S: Sensitivity

I: Electric current
from [electrometer 2](#)
(corresponds to high voltage)

C: Particle concentration
from [CPC](#)



UNIVERSITY OF MINNESOTA

2. Selection of appropriate voltage (spheres)



- **Low voltage: 20V (default voltage):** only removing excess ions; covering the whole size range
- **High voltage: Which voltage should we choose?**


I1	I2	Fitting equation	R ²
20V	2400	$S = 3.23 \times 10^7 I_2 - 5.4 \times 10^5 I_1$	0.9868
20V	1800	$S = 2.45 \times 10^7 I_2 - 1.49 \times 10^6 I_1$	0.9903
20V	1200	$S = 1.84 \times 10^7 I_2 - 2.5 \times 10^6 I_1$	0.9901
20V	600	$S = 1.89 \times 10^7 I_2 - 7.04 \times 10^6 I_1$	0.9661
20V	300	$S = 2.53 \times 10^7 I_2 - 1.55 \times 10^7 I_1$	0.9642

1200V is selected, which is a tradeoff between:

- ❑ **minimal sensitivity limit** (enough signal-to-noise ratio)
- ❑ **Goodness of fitting** (under different high voltages)



3. Validation (New method vs. SMPS) setup




- **High voltage: 1200V, Ion trap: 20V**
- **Compared to SMPS (Polydisperse particles)**



3. Validation A (New Method vs. SMPS)

□ Polydisperse KCl aerosol (salt)




Continuous Sampling



UNIVERSITY OF MINNESOTA

3. Validation B (New Method vs. SMPS)

□ Polydisperse **DEHS** aerosol (oil)





Continuous Sampling



UNIVERSITY OF MINNESOTA

4. Setup diagram- Agglomerate



- ❑ Furnace 1: **1200°C** (generating silver spheres); Furnace 2: **room temperature to 600°C** (sintering to tune aggregates structure)^a





^a Seong C. Kim, et al. (2009) Structural property effect of nanoparticle agglomerates on particle penetration through fibrous filter, *Aerosol Science and Technology*, 43:4, 344-355.



5. Results -Agglomerate

- Effect of sintering temperature (Change of particle morphology)




- Increasing temperature, particle size distribution changes.
- Increasing temperature, the sensitivity decreases significantly
- Comparing to spheres, higher sensitivities were found for aggregates




Outline

- ☐ Background and objective
- ☐ Methodology
- ☐ Experiments
- ☒ Summary





Summary



- ❑ The developed method (by combining one ion trap and one ESP) can successfully measure the Particle Geometric Surface Area of polydisperse aerosols with the following features:
 - With 1s time resolution
 - In a wide concentration range
 - Adapting to polydisperse particle distribution
 - Cost-effective
 - With simple setup



Future work


- Developing the method that can measure GSA of Agglomerate particles with 1s response time
- Adding APM (Aerosol Particle Mass Analyzer) to the set-up to offer more substantial and valuable information for agglomerates
- Investigating the effect of humidity on the GSA measurement by the developed method to apply it to different operating conditions




Thank you




Appendix 1-Practical model



Overall




Back side



Front side





Appendix-Calibration setup



- **Monodisperse Particle loss calibration:** The first NSAM charger and ion trap **on** while the second NSAM charger and ion trap **off**. Comparing results between two NSAMs.
- **Polydisperse calibration** is similar to Monodisperse calibration



Appendix 3-Loss calibration



$$\frac{\text{first current} - \text{second current}}{\text{first current}} < 5\%$$

The developed method has a negligible particle loss **in 40 to 300nm**, confirming the feasibility of new set-up

